JCI The Journal of Clinical Investigation

Studies of T- and B-Lymphocytes in Patients With Connective Tissue Diseases

J Clin Invest. 1973;52(7):3f-3f. https://doi.org/10.1172/JCI107184E1.

Erratum





Meet The Family

MICRO...FAST...DEPENDABLE...OBSOLESCENCE-FREE, MODULAR DESIGN instruments that pipette, dilute, dispense, read, count. And with greater precision and reproducibility than other comparable instruments.

All Micromedic Systems liquid delivery instruments offer: ■

- Repeatability of 0.5% C.V. or better.
- Continuous micrometer volume adjustment from 0 to 100% of pump capacity.

Automated Pipetting Station

- Eliminates manual effort in RIA.
- Offers exclusive Slurry Dispenser that automates separations using charcoal and other slurries.

MS 2 Spectrophotometer A

- Twin 17 ul cuvettes.
- Two seconds from feed to read.
- Stable: drift less than 0.001A per day.
- Repeatable: ± 0.001A or better.

- Eliminates tube handling by using standard MSI racks.
- 588-tube capacity.
- Multiple-user operation: each of 3 trays can be independently programmed ... even for isotope selection.
- Automatic shutdown.
- Interruption of automatic mode for manual counting without loss of index.

MICROMEDIC SYSTEMS, INC. Rohm and Haas Bldg., Phila., Pa. 19105 I would like salesman to call					
			Name	Date	
			Organization		
Address					
City	State	Zip			
Phone: Area Code	(Number)				

