Supplemental Data

Supplemental Figures




Figure S1. H&E stained femurs (H&E, 63x) are shown from 700 cGy-irradiated mice at day +10, treated every other day with saline or PTN.




Figure S2. Mean percentage p-ERK1/2 levels are shown in BM KSL cells from PTN $^{+/+}$ and PTN $^{-/-}$ mice (*P=0.005, n=3-5).

Supplemental Methods

Radiation Survival Studies and Hematopoietic Cell Transplantation Studies

Ten to 12 week old female C57Bl6 mice were irradiated with an LD50/30 radiation dose

(700 cGy) of TBI using a Cs137 irradiator. The first group of mice received intraperitoneal injection of 2 or 4 µg PTN or saline control in a volume of 100 ul starting at 24 hours post-irradiation. Treatment was administered every other day through day +14. For all other survival studies, mice received PTN subcutaneously (SQ). For the survival studies which compared PTN treatment versus PTN + Tipifarnib, PTN 4 µg SQ and Tipifarnib 1 mg (or vehicle), via gavage, were administered every other day through day +21. For hematopoietic cell transplant studies, 10-12 week old female C57Bl6 mice were irradiated with an LD100/30 dose of 850 cGy. Mice were then transplanted via tail vein with non-irradiated donor C57Bl6 BM cells (1 x 10⁵ cells/mouse). Mice then received PTN or saline beginning at 24 hours post-irradiation. Treatment was administered every other day through day +14.

HSC culture and Ras inhibition studies

BM HSCs were isolated and purified from mice as previously described (22). Lineage⁻ cells were stained with FITC-conjugated anti-CD34 (eBioscience), PE-conjugated anti-sca-1, and APC-conjugated anti-ckit (Becton Dickinson [BD]), or isotype controls. Sterile cell sorting was conducted on a BD FACS-Aria cytometer. Purified ckit⁺sca-1⁺lin⁻ (KSL) cells or CD34⁻c-kit⁺sca-1⁺lin⁻ (34⁻KSL) cells were collected into IMDM + 10% FBS + 1% pcn/strp.

In vitro inhibition of Ras signaling was performed by addition of the farnesyltransferase inhibitor FTI-276 (R&D Systems) at a concentration of 20 uM. In vivo inhibition of Ras signaling was accomplished by treating mice with the Ras inhibitor, Tipifarnib (SelleckChem) at +1 hr following TBI and every other day through day +21. Tipifarnib was dissolved in a 20% (2-Hydroxypropyl)-β-cyclodextrin solution at a concentration of 1 mg/100 μL. Mice were gavaged with 100 uL of Tipifarnib or cyclodextrin vehicle.

Competitive Repopulating Unit (CRU) Assays

BM 34 KSL cells from donor CD45.1 B6.SJL mice were isolated. The cells were cultured for 7 days in TSF media (IMDM containing 10%FBS, 1% pen-strep, 20 ng/ml Thrombopoietin, 125 ng/ml SCF, 50 ng/ml Flt-3 ligand, R&D) alone, TSF media + 100 ng/ml PTN, or TSF media + PTN + 20 uM FTI-276. Recipient CD45.2 C57Bl6 mice received 950 cGy total body irradiation (TBI) using a Cs137 irradiator and were then injected via tail vein with the progeny of 10 BM 34 KSL cells from the above culture conditions. Non-irradiated host BM MNCs (1 x 10⁵ cells/mouse) were injected as competitor cells. Multi-lineage hematopoietic engraftment was measured in the peripheral blood (PB) by flow cytometry, as previously described (16, 22).

Colony Forming Cell (CFC) Assays

Total CFC assays and colony-forming unit-granulocyte, erythroid, monocyte, megakaryocyte (CFU-GEMM) assays, to measure multipotent hematopoietic progenitor cells, were performed as we have previously described (8). For the studies of irradiated and transplanted mice, 2 x 10⁴ BM cells from each mouse were plated per CFC dish. For the in vitro studies, the progeny of 200 KSL cells cultured for 7 days were plated per CFC dish. Colonies were grown in MethoCult media (StemCell Technologies) for 14 days.

Flow Cytometric Analyses

KSL cells were cultured in the following conditions: IMDM media alone, IMDM + 100 ng/ml PTN (R&D), IMDM + 100 ng/ml PTN (R&D) + 20 uM FTI-276 (R&D). After 6 hours of treatment, cell cycle analysis was performed as previously described (8). After fixation and treatment with Saponin for permeabilization, cells were stained for 30 min with Ki-67 FITC antibody (BD), rinsed again, and then re-suspended in Saponin buffer with 7AAD (BD). At +24, +48, and +72

time points, Annexin V-FITC /7AAD staining was also performed to assess cell death per the manufacturers protocol (BD).

For phospho-flow cytometric analysis of pGrb2 (pS159, Sigma Aldrich SAB4504672), pERK1/2 (pT202/pY204, BD Biosciences 612952), and pALK (Y1604, Cell Signaling Technology 3341), the following methodology was used. KSL cells were stimulated with 100 ng/ml rPTN for 15-30 minutes at 37°C. Cells were fixed with BD Cytofix™ buffer (Catalogue number 554655) for 10 minutes at 37°C and then permeabilized by adding BD™ Phosflow Perm Buffer III (Catalogue number 558050) for 30 minutes on ice. Cells were washed twice and stained with the primary antibody for 1 hour. The cells were then rinsed and stained with the secondary antibodies prior to FACS analysis.

Immunofluorescence and Immunohistochemical Analyses

Hematoxylin and Eosin (H&E) staining was performed on femur sections from mice. Femurs were fixed overnight in 4% paraformaldehyde, decalcified and embedded in OCT medium. Cryostat sections were stained for H&E as previously described (8, 16). pERF (pT526, Abcam ab59178) analysis was performed by treating lineage negative BM cells for 30 min in IMDM or IMDM + 100 ng/ml rPTN. Cells were then spun onto glass slides using a Cytospin centrifuge at a density of 2 x 10⁴ cells/slide. Samples were fixed with 4% paraformaldehyde for 20 minutes, washed with PBS + 0.1% Tween, and blocked with PBS + 0.1% Tween + 20% goat serum. The primary antibody was diluted 1:100 in PBS + 0.1% Tween and added for 1 hour. After several rinses, the secondary antibody (goat anti-rabbit Alexa Fluor 488; Life Technologies) was added and the slides were rinsed and mounted with Vectashield Mounting Media with DAPI. The number of cells in each image which were pERF positive were counted.